Coussin thermique GRANDER®

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GRANDER GmbH.
Bergwerksweg 10, 6373 Jochberg, Austria
www.grander.com

GRANDER® sleeve

Experience cosy warmth or soothing cooling with the GRANDER® sleeve. Made of cotton and stretch fabric, it incorporates a gel pad containing the GRANDER® active medium. Its adjustable hook-and-loop fastener enables a customised fit for maximum wearing comfort, while its machine-washable cover prevents cold and heat loss.

Although the GRANDER® sleeve works at room temperature, it can also deliver a pleasant heating or cooling effect that last several hours – simply warm the gel compress inside the sleeve in hot water or place it in the fridge to cool. The eco-friendly gel can be used up to 400 times.

The sleeve can be used to treat a wide range of conditions such as tension, bruising, swelling and joint pain. It can be worn in various places including your neck, head, shoulder, arm, hip, haunches, pelvic area, legs and feet. The sleeve invigorates the body directly thanks to its integrated GRANDER® revitalisation effect and can help alleviate pain or tension.

To use for warming:

Heat in a pan:

Bring water to a boil and then take off the heat. Place the gel compress in the hot water for five to ten minutes. Carefully remove the hot compress from the water and insert it into the fabric cover. The GRANDER® sleeve gently moulds to your body and releases a pleasant, soothing warmth that relieves tension and promotes blood circulation – just perfect after a stressful day or while enjoying a relaxing moment of wellness.

Please note:

- > Your GRANDER® sleeve is designed to give off a pleasant level of warmth rather than intense heat. Do not overheat! (max. 50 °C/120 °F)
- > Do not heat the gel compress in the microwave!

To use for cooling:

Cool in the fridge or freezer:

The gel compress can be cooled before use, either for one to two hours in the freezer or for several hours in the fridge. It will then produce a pleasant cooling effect that guarantees relief, especially following sporting injuries or swelling or simply on a hot day.

Please note:

- > Do not use cold for more than 20 minutes!
- > Only use for longer if instructed to do so by your doctor.
- > Please consult your doctor or pharmacist if you suffer from vasoconstriction (narrowing of the blood vessels) or haemodynamic instability (irregular circulation).

Design: stretch fabric and cotton with length-adjustable hook-and-loop fastener

Machine-washable: at 30 °C (85 °F); do not tumble dry

Permitted temperature range: +50 °C (120 °F) down to –12 °C (10 °F)

N.B.:

- > Stop using the product if it is damaged.
- > Do not let babies or infants play with the gel compress.
- > Do not place a hot or cold gel compress directly on the skin.
- > There is a risk of scalding at very high temperatures and frostbite at very low temperatures.

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